



HEALTH X SCHOLARS



Health Scholars How to Use Quest 2

➤ Kit will Include:

- Quest 2 w/Controllers





Quest 2 Setup

- Unpack Quest 2 and set on table
- Power on Quest 2
 - Hold down power the button on the right side of the Quest 2.
- Grab 1 Controller
 - Only 1 controller is needed for setup.
- Put on Quest 2 Headset
 - Adjust headset straps as needed before putting on head.


<https://youtu.be/ZBSiZ5Pcjjg?t=175>





Quest 2 Room Boundary Setup

- Putting on the headset will start the room boundary Setup.
- A room boundary is needed for all Health Scholars Apps.
- The Quest 2 will always prompt for guardian setup in a new environment.
- The Quest 2 requires a 6ft x 6ft well-lit area away from windows and direct sunlight.






Room Boundary Setup

- Touch the floor softly with the controller in your hand.
- You will see a window asking you to continue to set the room boundaries.
- Pull the trigger on the controller to draw your play area.
- Make sure you have enough space between the play area and obstructions in your room. (walls, desks, chairs, etc.)
- Once you finish drawing the play area it will show you a grid.
- Observe this grid and walk around the room to make sure you have enough space.




How to Access Health Scholars Apps


- After the Guardian is created, the Quest 2 will go into the Health Scholars home environment.
- This is where a user will access the menu to access Health Scholars apps.
- The next step is to connect the Quest 2 to a Wi-Fi Network
 - The Quest 2 requires Wi-fi for use.






Quest 2 Wi-Fi Setup


- To connect the Quest 2 to Wi-Fi, use the controller to click the  icon.
- The Wi-Fi icon can be found on the right side of the menu.
- This will open the Wi-Fi menu.
- A window will appear showing available Wi-Fi networks.
- Use the controller to select the correct network and login.






Room Boundary cont.


- At any time the room boundary can be adjusted.
 - This helps for taller users.
- Use the controller to click on the  icon to change the room boundary.
 - A window will appear showing available settings.
- Use the controller to perform the Guardian setup.
- Both Stationary and Roomscale options are available for the Guardian Settings





OFB Quest 2 Screen Casting


- To cast what a user is seeing in VR, use the controller to select the  icon
- A window will appear showing available settings.
- Use the controller to select the correct device you want to cast to.
- Both devices will need to be on the same Wi-Fi network.





Health Scholars Apps

- The Health Scholars VR Apps can be found on the Quest 2's home screen.
- Use the controller to select the Health Scholars app of your choice.
- Health Scholars VR Apps will have a logo with the name of the program underneath.
- After starting the app click demo to start the VR simulation.





Quest 2 Power and Charging

- Keep the Quest 2 powered off and charging when not in use.
 - This prevents the headset from staying on because it thinks someone is wearing it.
- Hold down the power button until the Quest 2 powers off.
 - Hold down the power button for longer than 5 seconds

